



INSULATED MAILERS

Easy to Use | Sustainable | Dim Weight Reduction




Why Coldcorr?

- Material Cushioning Provides Superior Payload Protection Over Traditional Bubble Mailers Resulting in Less Damage During Transit
- Sizes Designed to Fit in Standard Carrier Overnight Outer Boxes
- Perfect for Smaller Items Like GLP-1 Medications to Help Reduce or Even Eliminate Dimensional Weight Charges
- Designed for Use With Both Refrigerated and Room Temperature Items
- Multiple Thickness Options Allow for Full Customization



Recycling Info

- Insulation Layer Manufactured From Post-Industrial Recycled Cotton and Other Natural Fibers
- Inner Component is Recyclable and Will Biodegrade in Time Under Normal Landfill Conditions
- Product is Approved Through the USDA Biopreferred Program
- White Poly Wrap and Outer Mailer Recyclable With #4 Plastics





COLD CORR®
ARIZONA | KENTUCKY | TENNESSEE

INSULATED MAILERS

Easy to Use | Sustainable | Dim Weight Reduction




Duration
12-36+ Hours

Test Profile
MOD. ISTA 7E

Temp Range
2-8 C

Coldcorr® Insulated Packaging Insulated Mailers

THICKNESS	OUTSIDE DIMS	INSIDE DIMS	QTY PER CASE	QTY PER PALLET
0.5"	9x12	7x10	150	1800
0.5"	12x15.5	10x13	100	1200
0.5"	14x19	12x17	70	840
1.0"	9x12	6x9	70	900
1.0"	12x15.5	9x12	48	576
1.0"	14x19	11x15	30	360
1.0"	19x24	16x20	20	240
1.5"	12x15.5	8x11	40	480
1.5"	15.5x12	11x8	40	480
1.5"	14x19	10x14	30	360
1.5"	19x14	14x10	30	360
1.5"	19x24	14x19	16	192



Testing Performed With 50ml Simulated Payload in a Chamber Environment